

Oil Rotary Vacuum Pumps Optional Suction and Exhaust Pipes

Suction Pipe

| Type | Products | Material | G-5DA | G-10DA G-20DA G-50DA | G-100D | GHD-030 | GHD-060 GHD-120 GLD-051 GLD-040 GCD-051X | GLD-136A GLD-136C GLD-201A GLD-201B GCD-136X GCD-201X | GLD-280A GLD-280B | Photo number |
|------------------|---|-------------|-------|----------------------------|--------|---------|--|--|----------------------|--------------|
| Hose port type | Hose port suction pipe (dia.8 x M20) | BS + plate | | | + 40 | | + 40 | + 50 | | |
| | Hose port suction pipe (dia.12 x M20) | BS + plate | | | + 40 | | + 40 | + 50 | | |
| | Hose port suction pipe (dia.15 x M20) | BS + plate | | | + 40 | | + 40 | + 50 | + 50 | |
| | Hose port suction pipe (dia.18 x M20) | BS + plate | | | + 40 | | + 40 | + 50 | + 50 | |
| | Hose port suction pipe (dia.22 x M20) | BS + plate | | | + 40 | | + 40 | + 50 | + 50 | |
| | Hose port suction pipe (dia.27 x 40) | ZDC + plate | | | | | | | | |
| | Hose port suction pipe (dia.27 x 50) | ZDC + plate | | | | | | | | |
| | Hose port suction pipe (dia.30 x 50) | ZDC + plate | | | | | | | | |
| For GHD-030 only | Hose port suction pipe for GHD-030 (dia.12) | BS + plate | | | | | | | | |
| | Hose port suction pipe for GHD-030 (dia.18) | BS + plate | | | | | | | | |
| | Suction pipe for GHD-030 (KF-16) | BS + plate | | | | | | | | |
| KF-25 type | KF-25 suction pipe (KF-25 x M20) | BS + plate | | | + 40 | | + 40 | + 50 | + 50 | |
| | KF-25 suction pipe (KF-25 x 40) | PPS | | | | | | | | |
| | KF-25 suction pipe (KF-25 x 50) | PPS | | | | | | | | |
| 1/4 screw type | A-type suction pipe (R1/4-L20 Female x M20) | BS + plate | | | + 40 | | + 40 | + 50 | + 50 | |
| | B-type suction pipe (R1/4 Male x M20) | BS + plate | | | + 40 | | + 40 | + 50 | + 50 | |
| | C-type suction pipe (R1/4-L10 Female x M20) | BS + plate | | | + 40 | | + 40 | + 50 | + 50 | |
| Adaptor | 40 adaptor (M20 Female) | BS + plate | | | | | | | | |
| | 50 adaptor (M20 Female) | BS + plate | | | | | | | | |

: Standard type : Replaceable (no adaptors necessary) + 40 : Replaceable (40 adaptor required) + 50 : Replaceable (50 adaptor required)

Suction Pipe



Exhaust Pipe



Exhaust Pipe

| Type | Products | Material | GLD-040 G-5DA G-10DA G-20DA G-50DA GLD-051 | G-100D | GHD-030 | GHD-060 GHD-120 GLD-136A GLD-136C GLD-201A GLD-201B GLD-280A | GCD-051X | GCD-136X GCD-201X | Photo number |
|----------------|--|------------|---|--------|---------|--|----------|----------------------|--------------|
| Hose port type | Hose port exhaust pipe (dia.12 x G3/4) | BS + plate | | | | | | | |
| | Hose port exhaust pipe (dia.15 x G3/4) | BS + plate | | | | | | | |
| | Hose port exhaust pipe (dia.18 x G3/4) | BS + plate | | | | | | | |
| | Hose port exhaust pipe (dia.27 x G1) | BS + plate | | | | | | | |
| KF flange type | KF-25 exhaust pipes (KF-25 x G3/4) | PPS | | | | | | | |
| | KF-25 exhaust pipes (KF-25 x G1) | PPS | | | | | | | |
| | KF-16 (G3/4) | BS + plate | | | | | | | |
| Exhaust pipe | Exhaust pipes Assy (G3/4) | PA | | | | | | | |
| | Exhaust pipes Assy (G1) | PA | | | | | | | |

Remarks: Different anti-corrosive quality as standard exhaust pipes

: Equipped with vacuum pump

: Replaceable (no adaptors necessary)